


## Solutions for Hand Fracture



□ Locking Compression Plate-small Y Shape (1.5)




**Screws:**  $\Phi$  1.5 Self- tapping locking screw- 3001815xx  
 $\Phi$  1.5 Cortex screw- 67115xx

**Indication:** It is used for the fracture on the top of phalanx.

Ref #	Size	Length	Material
08.311.01.20025802a	2Holes	17.5mm	A
08.311.01.20025803a	3Holes	22.5mm	
08.311.01.20025804a	4Holes	27.5mm	
08.311.01.20025805a	5Holes	32.5mm	
08.311.01.20025806a	6Holes	37.5mm	
08.311.01.20025807a	7Holes	42.5mm	
08.311.01.20025808a	8Holes	47.5mm	

□ Reconstruction Locking Compression Plate-1-small T Shape (1.5)



**Screws:**  $\Phi$  1.5 Self- tapping locking screw- 3001815xx  
 $\Phi$  1.5 Cortex screw- 67115xx

**Indication:** It is used for the fracture on the bottom of phalanx.

Ref #	Size	Length	Material
08.311.01.20035802a	2Holes	14.7mm	A
08.311.01.20035803a	3Holes	19.7mm	
08.311.01.20035804a	4Holes	24.7mm	
08.311.01.20035805a	5Holes	29.7mm	
08.311.01.20035806a	6Holes	34.7mm	
08.311.01.20035807a	7Holes	39.7mm	

### Mini Condylar Locking Compression Plate- 2



**Screws:**  $\Phi$ 2.0 Self- tapping locking screw- 3001820xx  
 $\Phi$ 2.0 Cortex screw- 67120xx

**Indication:** It is used for the fracture on metacarpal bone.

Ref #	Size	Length	Material
08.311.01.100970804a	4Holes	34mm	A
08.311.01.100970805a	5Holes	40mm	

### Medial/lateral Column Constructive Locking Plate(2.0)



**Screws:**  $\Phi$ 2.0 Self- tapping locking screw- 3001820xx  
 $\Phi$ 2.0 Cortex screw- 67120xx

**Indication:** For hand bone fixation.

Ref #	Size	Length	Material
08.311.01.100970802c	2Holes	19.5mm	A
08.311.01.100970803c	3Holes	26.3mm	
08.311.01.100970804c	4Holes	33.1mm	
08.311.01.100970805c	5Holes	39.9mm	
08.311.01.100970806c	6Holes	46.7mm	

Mini Locking Compression Plate-1



**Screws:**  $\Phi$ 2.0 Self- tapping locking screw- 3001820xx  
 $\Phi$ 2.0 Cortex screw- 67120xx

**Indication:** It is used for the fracture on metacarpal bone.

Ref #	Size	Length	Material
08.311.01.100080804a	4Holes	31mm	A
08.311.01.100080805a	5Holes	38mm	
08.311.01.100080806a	6Holes	45mm	
08.311.01.100080807a	7Holes	52mm	
08.311.01.100080808a	8Holes	59mm	

Reconstruction Locking Compression Plate-2-small T shape



**Screws:**  $\Phi$ 2.0 Self- tapping locking screw- 3001820xx  
 $\Phi$ 2.0 Cortex screw- 67120xx

**Indication:** For the base of the metacarpal bone.

Ref #	Size	Length	Material
08.311.01.200045802	2Holes	18.5mm	A
08.311.01.200045803	3Holes	25.3mm	
08.311.01.200045804	4Holes	32.1mm	
08.311.01.200045805	5Holes	38.9mm	
08.311.01.200045806	6Holes	45.7mm	
08.311.01.200045807	7Holes	52.5mm	



## Mini Locking Compression Plate-2 (2.0)



**Screws:**  $\Phi$ 2.0 Self- tapping locking screw- 3001820xx  
 $\Phi$ 2.0 Cortex screw- 67120xx

**Indication:** Used for fixation of the metacarpal bone or the foot AKIN procedure.

Ref #	Size	Length	Material
08.311.01.200015802	2Holes	11.8mm	A
08.311.01.200015803	3Holes	18.6mm	
08.311.01.200015804	4Holes	25.4mm	
08.311.01.200015805	5Holes	32.2mm	
08.311.01.200015806	6Holes	39mm	
08.311.01.200015807	7Holes	45.8mm	
08.311.01.200015808	8Holes	52.6mm	
08.311.01.200015809	9Holes	59.4mm	
08.311.01.200015810	10Holes	66.2mm	
08.311.01.200015811	11Holes	73mm	
08.311.01.200015812	12Holes	79.8mm	

## 1.5mm Self- tapping Locking Screw



**Matching Plate:** Used for 1.5mm locking plate

**With drill bit:**  $\Phi$  1.1 (WG REF # 30900700)

**Matching screwdriver:** wrench T5 (WG REF #31105107)

**Torque limited :** 0.25Nm (WG REF # 31103601)

Ref #	Size	Material
08.311.01.300181506	1.5×6mm	T
08.311.01.300181507	1.5×7mm	
08.311.01.300181508	1.5×8mm	
08.311.01.300181509	1.5×9mm	
08.311.01.300181510	1.5×10mm	
08.311.01.300181511	1.5×11mm	
08.311.01.300181512	1.5×12mm	
08.311.01.300181514	1.5×14mm	
08.311.01.300181516	1.5×16mm	
08.311.01.300181518	1.5×18mm	
08.311.01.300181520	1.5×20mm	
08.311.01.300181522	1.5×22mm	
08.311.01.300181524	1.5×24mm	

## Ø2.0mm Self- tapping Locking Screw



**Matching Plate:** It is used for 2.0mm locking plate.

**With drill bit:**  $\Phi$  1.5 (WG REF #30900800)

**Matching screwdriver:** wrench T6 (WG REF #31105106)

**Torque limited :** 0.4Nm (WG REF # 31103600)

Ref #	Size	Material
08.311.01.300182008	2.0×8mm	T
08.311.01.300182009	2.0×9mm	
08.311.01.300182010	2.0×10mm	
08.311.01.300182011	2.0×11mm	
08.311.01.300182012	2.0×12mm	
08.311.01.300182014	2.0×14mm	
08.311.01.300182016	2.0×16mm	
08.311.01.300182018	2.0×18mm	
08.311.01.300182020	2.0×20mm	
08.311.01.300182022	2.0×22mm	
08.311.01.300182024	2.0×24mm	
08.311.01.300182026	2.0×26mm	
08.311.01.300182028	2.0×28mm	
08.311.01.300182030	2.0×30mm	
08.311.01.300182032	2.0×32mm	

## Ø1.5mm Cortex Screw



**Matching Plate:** It is used for 2.0mm locking plate.


**With drill bit:**  $\Phi$  1.5 (WG REF #30900700)

**Matching screwdriver:** wrench SW 1.5 (WG REF #31105101)

Ref #	Size	Material
08.301.11.6711505	1.5×5mm	T
08.301.11.6711506	1.5×6mm	
08.301.11.6711507	1.5×7mm	
08.301.11.6711508	1.5×8mm	
08.301.11.6711509	1.5×9mm	
08.301.11.6711510	1.5×10mm	
08.301.11.6711511	1.5×11mm	
08.301.11.6711512	1.5×12mm	
08.301.11.6711514	1.5×14mm	
08.301.11.6711516	1.5×16mm	
08.301.11.6711518	1.5×18mm	
08.301.11.6711520	1.5×20mm	
08.301.11.6711522	1.5×22mm	
08.301.11.6711524	1.5×24mm	
08.301.11.6711526	1.5×26mm	
08.301.11.6711528	1.5×28mm	
08.301.11.6711530	1.5×30mm	
08.301.11.6711532	1.5×32mm	
08.301.11.6711534	1.5×34mm	
08.301.11.6711536	1.5×36mm	
08.301.11.6711538	1.5×38mm	
08.301.11.6711540	1.5×40mm	



## Ø2.0mm Cortex Screw



**Matching Plate:** It is used for 2.0mm locking plate.

**With drill bit:**  $\phi$  1.5 (WG REF #30900800)

**Matching screwdriver:** wrench SW1.5 (WG REF #31105102)

Ref #	Size	Material
08.301.11.6712005	2.0×5mm	T
08.301.11.6712006	2.0×6mm	
08.301.11.6712007	2.0×7mm	
08.301.11.6712008	2.0×8mm	
08.301.11.6712009	2.0×9mm	
08.301.11.6712010	2.0×10mm	
08.301.11.6712011	2.0×11mm	
08.301.11.6712012	2.0×12mm	
08.301.11.6712014	2.0×14mm	
08.301.11.6712016	2.0×16mm	
08.301.11.6712018	2.0×18mm	
08.301.11.6712020	2.0×20mm	
08.301.11.6712022	2.0×22mm	
08.301.11.6712024	2.0×24mm	
08.301.11.6712026	2.0×26mm	
08.301.11.6712028	2.0×28mm	
08.301.11.6712030	2.0×30mm	
08.301.11.6712032	2.0×32mm	
08.301.11.6712034	2.0×34mm	
08.301.11.6712036	2.0×36mm	
08.301.11.6712038	2.0×38mm	
08.301.11.6712040	2.0×40mm	








www.wegortho.com

## Solutions for Foot and Ankle Fracture



## ☐ Astragalus Locking Plate




**Screws:**  $\Phi$ 2.4 Self- tapping locking screw- 3001924xx  
 $\Phi$ 2.7 Self- tapping locking screw- 3001927xx  
 $\Phi$ 2.4 Cortex screw- 61124xx

**Indication:** It is used for Astragalus fracture.

Ref #	Size	Length	Material
08.311.01.100950800	6Holes	38mm	A

## ☐ Scaphoid Locking Plate



**Screws:**  $\Phi$ 2.4 Self- tapping locking screw- 3001824xx  
 $\Phi$ 2.7 Self- tapping locking screw- 3001827xx  
 $\Phi$ 2.4 Cortex screw- 61124xx

**Indication:** It is used for Foot Scaphoid fracture.

Ref #	Size	Length	Material
08.311.01.100955800	Fixed Size	51mm	A

## ☐ Cuboid Locking Plate




**Screws:**  $\Phi$ 2.4 Self- tapping locking screw- 3001824xx  
 $\Phi$ 2.7 Self- tapping locking screw- 3001827xx  
 $\Phi$ 2.4 Cortex screw- 61124xx

**Indication:** It is used for cuboid fracture.

Ref #	Size	Length	Material
08.311.01.200051800	Left	26mm	A
08.311.01.200052800	Right	26mm	



### Reconstruction Locking Compression Plate-2- t Shape




**Screws:**  $\Phi$ 2.4 Self- tapping locking screw- 3001824xx  
 $\Phi$ 2.7 Self- tapping locking screw- 3001827xx  
 $\Phi$ 2.4 Cortex screw- 61124xx

**Indication:** It is used for metacarpal and metatarsal fracture.

Ref #	Size	Length	Material
08.311.01.200040802	2Holes	21.25mm	A
08.311.01.200040803	3Holes	28.5mm	
08.311.01.200040804	4Holes	35.75mm	
08.311.01.200040805	5Holes	43mm	
08.311.01.200040806	6Holes	50.25mm	
08.311.01.200040807	7Holes	57.5mm	

### Reconstruction Locking Compression Plate-1- T Shape (2.4)



**Screws:**  $\Phi$ 2.4 Self- tapping locking screw- 3001824xx  
 $\Phi$ 2.7 Self- tapping locking screw- 3001827xx  
 $\Phi$ 2.4 Cortex screw- 61124xx

**Indication:** It is used for metacarpal and metatarsal fracture.

Ref #	Size	Length	Material
08.311.01.200030802	2Holes	21.55mm	A
08.311.01.200030803	3Holes	28.8mm	
08.311.01.200030804	4Holes	36.05mm	
08.311.01.200030805	5Holes	43.3mm	
08.311.01.200030806	6Holes	50.55mm	
08.311.01.200030807	7Holes	57.8mm	

## Reconstruction Locking Compression Plate-Y Shape




**Screws:**  $\Phi$ 2.4 Self- tapping locking screw- 3001824xx  
 $\Phi$ 2.7 Self- tapping locking screw- 3001827xx  
 $\Phi$ 2.4 Cortex screw- 61124xx

**Indication:** It is used for metacarpal and metatarsal fracture.

Ref #	Size	Length	Material
08.311.01.200020802	2Holes	23.75mm	A
08.311.01.200020803	3Holes	31mm	
08.311.01.200020804	4Holes	38.25mm	
08.311.01.200020805	5Holes	45.5mm	
08.311.01.200020806	6Holes	52.75mm	
08.311.01.200020807	7Holes	60mm	

## 2.4mm Condylar Locking Compression Plate



**Screws:**  $\Phi$ 2.4 Self- tapping locking screw- 3001824xx  
 $\Phi$ 2.7 Self- tapping locking screw- 3001827xx  
 $\Phi$ 2.4 Cortex screw- 61124xx

**Indication:** For osteophyte fixation.

Ref #	Size	Length	Material
08.311.01.100975802	2Holes	22.45mm	A
08.311.01.100975803	3Holes	29.7mm	
08.311.01.100975804	4Holes	36.95mm	
08.311.01.100975805	5Holes	44.2mm	
08.311.01.100975806	6Holes	51.45mm	
08.311.01.100975807	7Holes	58.7mm	



### Mini Locking Compression Plate-1 (2.7)



**Screws:**  $\Phi$ 2.4 Self- tapping locking screw- 3001824xx  
 $\Phi$ 2.7 Self- tapping locking screw- 3001827xx  
 $\Phi$ 2.7 Cortex screw- 61127xx

**Indication:** It is used for metatarsal fracture.

Ref #	Size	Length	Material
08.311.01.100080804	4Holes	40mm	A
08.311.01.100080805	5Holes	49mm	
08.311.01.100080806	6Holes	58mm	
08.311.01.100080807	7Holes	67mm	
08.311.01.100080808	8Holes	76mm	
08.311.01.100080810	10Holes	94mm	

### Miniature Osteophyte Compression Locking Plate(2.7)




**Screws:**  $\Phi$ 2.4 Self- tapping locking screw- 3001824xx  
 $\Phi$ 2.7 Self- tapping locking screw- 3001827xx  
 $\Phi$ 2.4 Cortex screw- 61124xx

**Indication:** For osteophyte fixation.

Ref #	Size	Length	Material
08.311.01.100940802	2Holes	22.45mm	A
08.311.01.100940803	3Holes	29.7mm	
08.311.01.100940804	4Holes	36.95mm	
08.311.01.100940805	5Holes	44.2mm	
08.311.01.100940806	6Holes	51.45mm	
08.311.01.100940807	7Holes	58.7mm	

## Locking Compression Plate-small T Shape(2.7)



**Screws:**  $\Phi$ 2.4 Self- tapping locking screw- 3001824xx  
 $\Phi$ 2.7 Self- tapping locking screw- 3001827xx  
 $\Phi$ 2.7 Cortex screw- 61127xx

**Indication:** It is used for metatarsal fracture.

Ref #	Size	Length	Material
08.311.01.100965803	3Holes	32mm	A
08.311.01.100965804	4Holes	39.5mm	

## Locking Compression Plate-small L Shape(2.7)



**Screws:**  $\Phi$ 2.4 Self- tapping locking screw- 3001824xx  
 $\Phi$ 2.7 Self- tapping locking screw- 3001827xx  
 $\Phi$ 2.7 Cortex screw- 61127xx

**Indication:** It is used for metatarsal fracture.

Ref #	Size	Length	Material
08.311.01.100961803	3 holes/left	32.3mm	A
08.311.01.100962803	3 holes/right	32.3mm	
08.311.01.100961804	4 holes/left	39.8mm	
08.311.01.100962804	4 holes/right	39.8mm	

## Locking Compression Plate-small Oblique L Shape (2.7)




**Screws:**  $\Phi$ 2.4 Self- tapping locking screw- 3001824xx  
 $\Phi$ 2.7 Self- tapping locking screw- 3001827xx  
 $\Phi$ 2.7 Cortex screw- 61127xx

**Indication:** It is used for metatarsal fracture.

Ref #	Size	Length	Material
08.311.01.100969803	3 holes/left	34.4mm	A
08.311.01.100968803	3 holes/right	34.4mm	
08.311.01.100969804	4 holes/left	41.9mm	
08.311.01.100968804	4 holes/right	41.9mm	

## 1/3 Tubular Plate



**Screws:**  $\Phi$  3.5 Cortex screw- 61335xx

**Indication:** It is used for fibula fracture, olecranon, pelvis, scapula, clavicle.

Ref #	Size	Material
08.301.01.5048004	4Holes	T
08.301.01.5048005	5Holes	
08.301.01.5048006	6Holes	
08.301.01.5048007	7Holes	
08.301.01.5048008	8Holes	



### Ø2.4mm Self- tapping Locking Screw



- Matching Plate: It is used for 2.4/2.7mm locking plate.

---

- With drill bit:  $\phi$  1.8 (WG REF #30901602)

---

- Matching screwdriver: wrench T8 (WG REF #31105105)

---

- Torque limited : 0.8Nm (WG REF # 31101701)

### Ø2.7mm Self- tapping Locking Screw



- Matching Plate: It is used for 2.4/2.7mm locking plate.

---

- With drill bit:  $\phi$  2.0 (WG REF #30901601)

---

- Matching screwdriver: wrench T8 (WG REF #31105105)

---


- Torque limited : 0.8Nm (WG REF # 31101701)

Ref #	Size	Material
08.311.01.300182407	2.4×7mm	T
08.311.01.300182408	2.4×8mm	
08.311.01.300182409	2.4×9mm	
08.311.01.300182410	2.4×10mm	
08.311.01.300182411	2.4×11mm	
08.311.01.300182412	2.4×12mm	
08.311.01.300182414	2.4×14mm	
08.311.01.300182416	2.4×16mm	
08.311.01.300182418	2.4×18mm	
08.311.01.300182420	2.4×20mm	
08.311.01.300182422	2.4×22mm	
08.311.01.300182424	2.4×24mm	
08.311.01.300182426	2.4×26mm	
08.311.01.300182428	2.4×28mm	
08.311.01.300182430	2.4×30mm	
08.311.01.300182432	2.4×32mm	

Ref #	Size	Material
08.311.01.300182708	2.7×8mm	T
08.311.01.300182710	2.7×10mm	
08.311.01.300182712	2.7×12mm	
08.311.01.300182714	2.7×14mm	
08.311.01.300182715	2.7×15mm	
08.311.01.300182716	2.7×16mm	
08.311.01.300182717	2.7×17mm	
08.311.01.300182718	2.7×18mm	
08.311.01.300182719	2.7×19mm	
08.311.01.300182720	2.7×20mm	
08.311.01.300182722	2.7×22mm	
08.311.01.300182724	2.7×24mm	
08.311.01.300182726	2.7×26mm	
08.311.01.300182728	2.7×28mm	
08.311.01.300182730	2.7×30mm	
08.311.01.300182732	2.7×32mm	
08.311.01.300182734	2.7×34mm	
08.311.01.300182736	2.7×36mm	
08.311.01.300182738	2.7×38mm	
08.311.01.300182740	2.7×40mm	
08.311.01.300182742	2.7×42mm	
08.311.01.300182744	2.7×44mm	
08.311.01.300182746	2.7×46mm	
08.311.01.300182748	2.7×48mm	
08.311.01.300182750	2.7×50mm	
08.311.01.300182755	2.7×55mm	
08.311.01.300182760	2.7×60mm	



## Ø2.4mm Cortex Screw



**Matching Plate:** For 2.4mm locking plate.


**With drill bit:**  $\phi$  1.8 (WG REF #30901602)

**Matching screwdriver:** wrench SW2.0 (WG REF #31105100)

Ref #	Size	Material
08.301.11.6112405	2.4×5mm	T
08.301.11.6112406	2.4×6mm	
08.301.11.6112407	2.4×7mm	
08.301.11.6112408	2.4×8mm	
08.301.11.6112409	2.4×9mm	
08.301.11.6112410	2.4×10mm	
08.301.11.6112411	2.4×11mm	
08.301.11.6112412	2.4×12mm	
08.301.11.6112413	2.4×13mm	
08.301.11.6112414	2.4×14mm	
08.301.11.6112416	2.4×16mm	
08.301.11.6112418	2.4×18mm	
08.301.11.6112420	2.4×20mm	
08.301.11.6112422	2.4×22mm	
08.301.11.6112424	2.4×24mm	
08.301.11.6112426	2.4×26mm	
08.301.11.6112428	2.4×28mm	
08.301.11.6112430	2.4×30mm	
08.301.11.6112432	2.4×32mm	
08.301.11.6112434	2.4×34mm	
08.301.11.6112436	2.4×36mm	
08.301.11.6112438	2.4×38mm	
08.301.11.6112440	2.4×40mm	
08.301.11.6112442	2.4×42mm	
08.301.11.6112444	2.4×44mm	
08.301.11.6112446	2.4×46mm	
08.301.11.6112448	2.4×48mm	
08.301.11.6112450	2.4×50mm	



## Ø2.7mm Cortex Screw



Matching Plate: For 2.7mm locking plate.

With drill bit:  $\Phi$ 2.0 (WG REF #30901601)

Matching screwdriver: wrenchSW2.5(WG REF #31101302)

Ref #	Size	Material
08.301.11.6112706	2.7×6mm	T
08.301.11.6112708	2.7×8mm	
08.301.11.6112710	2.7×10mm	
08.301.11.6112712	2.7×12mm	
08.301.11.6112714	2.7×14mm	
08.301.11.6112715	2.7×15mm	
08.301.11.6112716	2.7×16mm	
08.301.11.611271	2.7×17mm	
08.301.11.6112718	2.7×18mm	
08.301.11.6112719	2.7×19mm	
08.301.11.6112720	2.7×20mm	
08.301.11.6112722	2.7×22mm	
08.301.11.6112724	2.7×24mm	
08.301.11.6112726	2.7×26mm	
08.301.11.6112728	2.7×28mm	
08.301.11.6112730	2.7×30mm	
08.301.11.6112732	2.7×32mm	
08.301.11.6112734	2.7×34mm	
08.301.11.6112736	2.7×36mm	
08.301.11.6112738	2.7×38mm	
08.301.11.6112740	2.7×40mm	
08.301.11.6112742	2.7×42mm	
08.301.11.6112744	2.7×44mm	
08.301.11.6112746	2.7×46mm	
08.301.11.6112748	2.7×48mm	
08.301.11.6112750	2.7×50mm	
08.301.11.6112752	2.7×52mm	
08.301.11.6112754	2.7×54mm	
08.301.11.6112756	2.7×56mm	
08.301.11.6112758	2.7×58mm	
08.301.11.6112760	2.7×60mm	



### Ø2.4mm Star Head Cortex Screw



**Matching Plate:** It is used for 2.4/2.7mm locking plate.


**With drill bit:**  $\phi$  1.8 (WG REF #08.301.22.31101602)

**Matching screwdriver:** T8 Star Screwdriver Shaft (WG REF #31105402)

**Torque limited :** 0.8Nm (WG REF # 31101701)

Ref #	Size	Material
08.301.13.6132410	2.4×10mm	T
08.301.13.6132412	2.4×12mm	
08.301.13.6132414	2.4×14mm	
08.301.13.6132416	2.4×16mm	
08.301.13.6132418	2.4×18mm	
08.301.13.6132420	2.4×20mm	
08.301.13.6132422	2.4×22mm	
08.301.13.6132424	2.4×24mm	
08.301.13.6132426	2.4×26mm	

### Ø2.7mm Star Head Cortex Screw



**Matching Plate:** It is used for 2.4/2.7mm locking plate.

**With drill bit:**  $\phi$  2.0 (WG REF #08.301.22.31101601)

**Matching screwdriver:** T8 Star Screwdriver Shaft (WG REF #31105402)

**Torque limited :** 0.8Nm (WG REF # 31101701)

Ref #	Size	Material
08.301.13.6132710	2.7×10mm	T
08.301.13.6132712	2.7×12mm	
08.301.13.6132714	2.7×14mm	
08.301.13.6132716	2.7×16mm	
08.301.13.6132718	2.7×18mm	
08.301.13.6132720	2.7×20mm	
08.301.13.6132722	2.7×22mm	
08.301.13.6132724	2.7×24mm	
08.301.13.61327126	2.7×26mm	