


Case Information:

Date: \_\_\_\_\_

Surgeon: \_\_\_\_\_


Patient ID / Booking No.: \_\_\_\_\_

**TOP TRAY – INSTRUMENTS**




(Theatre Usage Code)	Supplier Code	Number & Description	Assembly	Standard Qty	Case Usage
4362	30902500	1. Periosteal Elevator		1	
4361	30902400	2. Reduction Forceps With Points		2	
4400	31103000	3. Plate bender		2	
4386	31101701	4. 0.25Nm Torque Limiting Attachment (for 1.5mm locking screws)		1	
4405	31103600	5. 0.4Nm Torque Limiting Attachment (for 2.0mm locking screws)		1	
4368	30905900	6. Hohmann Retractors		1	
4462	31109400	Crab Claw Reduction Forceps (loose)		2	

**BOTTOM TRAY - INSTRUMENTS**




















Product Code (Theatre Usage Code)	Supplier Code	Number & Description	Assembly	Standard Qty	Case Usage
4365	30902900	7. Screwdriver Handle with Quick Coupling		1	
4401	31103400	8. Depth Gauge 25mm (Small)	4 PARTS	1	
4417	31105105	9. T6 Stardrive Screwdriver		1	
4418	31105106	10. T5 Stardrive Screwdriver		1	
	Blank	11. 1.5mm Hex Screwdriver Shaft (2.0mm Cortex)		1	
4364	30902700	12. Senn Miller Retract Hook (Cats Paws)		2	


2.0mm MODULE – SCREWS









Barcode	Supplier Code	Number & Description	Standard Qty	Case Usage
	08.301.11.6712008	1. Cortex screw 2.0x8mmTAV	4	
	08.301.11.6712009	1. Cortex screw 2.0x9mmTAV	4	
	08.301.11.6712010	1. Cortex screw 2.0x10mmTAV	4	
	08.301.11.6712011	1. Cortex screw 2.0x11mmTAV	4	
	08.301.11.6712012	1. Cortex screw 2.0x12mmTAV	4	
	08.301.11.6712014	1. Cortex screw 2.0x14mmTAV	4	
	08.301.11.6712016	1. Cortex screw 2.0x16mmTAV	4	
	08.301.11.6712018	1. Cortex screw 2.0x18mmTAV	4	
	08.311.01.300182008	2. Self-tapping Locking screw 2.0x8mm TAV	4	
	08.311.01.300182009	2. Self-tapping Locking screw 2.0x9mm TAV	4	
	08.311.01.300182010	2. Self-tapping Locking screw 2.0x10mm TAV	4	

 0106912547058540	08.311.01.300182011	2. Self-tapping Locking screw 2.0×11mm TAV	4	
 0106912547058557	08.311.01.300182012	2. Self-tapping Locking screw 2.0×12mm TAV	4	
 0106912547058564	08.311.01.300182014	2. Self-tapping Locking screw 2.0×14mm TAV	4	
 0106912547058571	08.311.01.300182016	2. Self-tapping Locking screw 2.0×16mm TAV	4	
 0106912547058588	08.311.01.300182018	2. Self-tapping Locking screw 2.0×18mm TAV	4	
 0106912547082613	08.301.11.6712020	3. Cortex screw 2.0x20mmTAV	4	
 0106912547082620	08.301.11.6712022	3. Cortex screw 2.0x22mmTAV	4	
 0106912547082637	08.301.11.6712024	3. Cortex screw 2.0x24mmTAV	4	
 0106912547082644	08.301.11.6712026	3. Cortex screw 2.0x26mmTAV	4	
 0106912547082651	08.301.11.6712028	3. Cortex screw 2.0x28mmTAV	4	
 0106912547082668	08.301.11.6712030	3. Cortex screw 2.0x30mmTAV	4	
 0106912547058595	08.311.01.300182020	3. Self-tapping Locking screw 2.0×20mm TAV	4	
 0106912547058601	08.311.01.300182022	3. Self-tapping Locking screw 2.0×22mm TAV	4	
 0106912547058618	08.311.01.300182024	3. Self-tapping Locking screw 2.0×24mm TAV	4	
 0106912547058625	08.311.01.300182026	3. Self-tapping Locking screw 2.0×26mm TAV	4	
 0106912547058632	08.311.01.300182028	3. Self-tapping Locking screw 2.0×28mm TAV	4	
 0106912547058649	08.311.01.300182030	3. Self-tapping Locking screw 2.0×30mm TAV	4	

**2.0mm MODULE – PLATES & INSTRUMENTS**













Barcode	Supplier Code	Number & Description	Standard Qty	Case Usage
 0106912547112945	08.311.01.100970805a	4. 2.0 mm Mini Condylar Locking Compression Plate-2 5 Holes Ti	1	
 0106912547113775	08.311.01.100080804a	5. 2.0 mm Locking Compression Plate 4 Holes Ti	2	
 0106912547113782	08.311.01.100080805a	5. 2.0 mm Locking Compression Plate 5 Holes Ti	2	
 0106912547113799	08.311.01.100080806a	5. 2.0 mm Locking Compression Plate 6 Holes Ti	2	
 0106912547113805	08.311.01.100080807a	5. 2.0 mm Locking Compression Plate 7 Holes Ti	1	
 0106912547113812	08.311.01.100080808a	5. 2.0 mm Locking Compression Plate 8 Holes Ti	2	
 0106912547056027	08.311.01.200035807	6. 2.0 mm Locking Reconstruction T-Plate (2 Hole) 7 Holes Ti	0	
 0106912547055181	08.311.01.200025807	7. 2.0 mm Locking Reconstruction Y-Plate 7 Holes Ti	0	
 0106912547056829	08.311.01.200045803	8. 2.0 mm Locking Reconstruction T-Plate (3 Hole) 3 Holes Ti	1	
 0106912547056836	08.311.01.200045804	8. 2.0 mm Locking Reconstruction T-Plate (3 Hole) 4 Holes Ti	2	
 0106912547056843	08.311.01.200045805	8. 2.0 mm Locking Reconstruction T-Plate (3 Hole) 5 Holes Ti	1	
 0106912547056850	08.311.01.200045806	8. 2.0 mm Locking Reconstruction T-Plate (3 Hole) 6 Holes Ti	1	

 0106912547056867	08.311.01.200045807	8. 2.0 mm Locking Reconstruction T-Plate (3 Hole) 7 Holes Ti	1	
 0106912547110828	08.311.01.100970803c	9. 2.0 mm Locking Condylar Plate 3 Holes Ti	1	
 0106912547110835	08.311.01.100970804c	9. 2.0 mm Locking Condylar Plate 4 Holes Ti	1	
 0106912547110842	08.311.01.100970805c	9. 2.0 mm Locking Condylar Plate 5 Holes Ti	1	
 0106912547110859	08.311.01.100970806c	9. 2.0 mm Locking Condylar Plate 6 Holes Ti	1	
	30901200	10. Drill Guide 2.0-1.5mm	1	
	Blank	12. 1.5mm Hex Screwdriver Shaft (2mm Cortex)	1	
	31107600	13. T6 Stardrive Screwdriver Shaft	2	
	30900502	14. Screw Tap for Cortex 2.0mm	1	
	30900102	14. Countersink for 1.5/2.0mm screws	1	
	31104200	15. 2.0mm Plate Bender	2	
	31103501	16. 1.5mm Screw in Drill Guide for 2.0mm Locking Screws	5	
<b>Barcode</b>	<b>Supplier Code</b>	<b>DISPOSABLE INSTRUMENTS (SINGLE USE ONLY)</b>	<b>Standard Qty</b>	<b>Case Usage</b>
 0106912547330015	30900800	11. 1.5mm Drill Bit	2	

1.5mm MODULE – SCREWS




Barcode	Supplier Code	Number & Description	Standard Qty	Case Usage
 0106912547082040	08.301.11.6711506	1. Cortex screw 1.5x6mm	6	
 0106912547082057	08.301.11.6711507	1. Cortex screw 1.5x7mm	6	
 0106912547082064	08.301.11.6711508	1. Cortex screw 1.5x8mm	6	
 0106912547082071	08.301.11.6711509	1. Cortex screw 1.5x9mm	6	
 0106912547082088	08.301.11.6711510	1. Cortex screw 1.5x10mm	6	
 0106912547082095	08.301.11.6711511	1. Cortex screw 1.5x11mm	6	
 0106912547082101	08.301.11.6711512	1. Cortex screw 1.5x12mm	6	
 0106912547082125	08.301.11.6711514	1. Cortex screw 1.5x14mm	6	









 0106912547082132	08.301.11.6711516	1. Cortex screw 1.5x16mm	6	
 0106912547082149	08.301.11.6711518	1. Cortex screw 1.5x18mm	6	
 0106912547120568	08.311.01.300181508	2. Self-tapping Locking screw 1.5x8mm	6	
 0106912547120575	08.311.01.300181509	2. Self-tapping Locking screw 1.5x9mm	6	
 0106912547120582	08.311.01.300181510	2. Self-tapping Locking screw 1.5x10mm	6	
 0106912547120599	08.311.01.300181511	2. Self-tapping Locking screw 1.5x11mm	6	
 0106912547120605	08.311.01.300181512	2. Self-tapping Locking screw 1.5x12mm	6	
 0106912547120612	08.311.01.300181514	2. Self-tapping Locking screw 1.5x14mm	6	
 0106912547120629	08.311.01.300181516	2. Self-tapping Locking screw 1.5x16mm	6	
 0106912547120636	08.311.01.300181518	2. Self-tapping Locking screw 1.5x18mm	6	



**1.5mm MODULE – PLATES & INSTRUMENTS**



Barcode	Supplier Code	Number & Description	Standard Qty	Case Usage
	08.311.01.200010804a	3. 1.5 mm Locking Compression Plates Straight 4 Holes	1	
n/a	08.311.01.200010805a	3. 1.5 mm Locking Compression Plates Straight 5 Holes	0	
	08.311.01.200010806a	3. 1.5 mm Locking Compression Plates Straight 6 Holes	1	
	08.311.01.200015802a	4. 1.5 mm Locking Plates Recon 2 Holes	1	
	08.311.01.200015803a	4. 1.5 mm Locking Plates Recon 3 Holes	1	
	08.311.01.200015804a	4. 1.5 mm Locking Plates Recon 4 Holes	1	
	08.311.01.200015805a	4. 1.5 mm Locking Plates Recon 5 Holes	1	
	08.311.01.200015806a	4. 1.5 mm Locking Plates Recon 6 Holes	1	
	08.311.01.200015807a	4. 1.5 mm Locking Plates Recon 7 Holes	1	

 0106912547114659	08.311.01.200015808a	4. 1.5 mm Locking Plates Recon 8 Holes	0	
 0106912547120261	08.311.01.20035802a	5. 1.5 mm Reconstruction Locking Compression T-Plates 2 Holes	1	
 0106912547120278	08.311.01.20035803a	5. 1.5 mm Reconstruction Locking Compression T-Plates 3 Holes	1	
 0106912547120285	08.311.01.20035804a	5. 1.5 mm Reconstruction Locking Compression T-Plates 4 Holes	1	
 0106912547120292	08.311.01.20035805a	5. 1.5 mm Reconstruction Locking Compression T-Plates 5 Holes	1	
 0106912547120308	08.311.01.20035806a	5. 1.5 mm Reconstruction Locking Compression T-Plates 6 Holes	1	
 0106912547120315	08.311.01.20035807a	5. 1.5 mm Reconstruction Locking Compression T-Plates 7 Holes	1	
n/a	30900102	6. Countersink for 1.5 screws	1	
n/a	31107801	7. 1.5mm Screw in Drill Guide	5	
n/a	30901100	8. Drill Guide 1.1-1.5mm	1	
n/a	31107900	9. 1.5mm Plate Bender	2	
n/a	31107700	10. T5 Stardrive Screwdriver Shaft – Quick Connect	2	
n/a	30900400	11. Screw Tap for Cortex 1.5mm	2	
<b>Barcode</b>	<b>Supplier Code</b>	<b>DISPOSABLE INSTRUMENTS (SINGLE USE ONLY)</b>	<b>Standard Qty</b>	<b>Case Usage</b>
 0106912547330011	30900700	12. 1.1mm Drill Bit	2	