



## Solutions for Clavicle Fracture



## □ Clavicle Locking Compression Plate-2



**Screws:**  $\Phi$ 3.5 Self- tapping locking screw- 1009835xx  
 $\Phi$ 3.5 Cortex screw- 61335xx  
 $\Phi$ 4.0 Cancellous screw- 61540xx

**Indication:** It is used for Clavicle Shaft fracture.

Ref #	Size	Length	Material
08.311.01.100034806	6 holes/left	94mm	A
08.311.01.100034807	7 holes/left	105mm	
08.311.01.100034808	8 holes/left	122mm	
08.311.01.100045806	6 holes/right	94mm	
08.311.01.100045807	7 holes/right	105mm	
08.311.01.100045808	8 holes/right	122mm	

## □ Clavicle Locking Compression Plate-2



**Screws:**  $\Phi$ 3.5 Self- tapping locking screw- 1009835xx  
 $\Phi$ 3.5 Cortex screw- 61335xx  
 $\Phi$ 4.0 Cancellous screw- 61540xx

**Indication:** Used for mid-column and lateral-column Clavicle Shaft fractures.

Ref #	Size	Length	Material
08.311.01.100034706	6 holes/left	94mm	T
08.311.01.100034707	7 holes/left	105mm	
08.311.01.100034708	8 holes/left	122mm	
08.311.01.100034709	9 holes/left	136mm	
08.311.01.100034710	10 holes/left	150mm	
08.311.01.100034712	12 holes/left	164mm	
08.311.01.100045706	6 holes/right	94mm	
08.311.01.100045707	7 holes/right	105mm	
08.311.01.100045708	8 holes/right	122mm	
08.311.01.100045709	9 holes/right	136mm	
08.311.01.100045710	10 holes/right	150mm	
08.311.01.100045712	12 holes/right	164mm	

### Clavicle Locking Compression Plate-3




**Screws:**  $\Phi$ 2.4 Self- tapping locking screw- 3001824xx  
 $\Phi$ 2.7 Self- tapping locking screw- 3001827xx  
 $\Phi$ 3.5 Self- tapping locking screw- 1009835xx  
 $\Phi$ 3.5 Cortex screw- 61335xx  
 $\Phi$ 4.0 Cancellous screw- 61540xx

**Indication:** It is used for fracture on the apex of clavicle or combined Shaft fracture.



Ref #	Size	Length	Material
08.311.01.100033803	3 holes/left	69mm	A
08.311.01.100033804	4 holes/left	81mm	
08.311.01.100033805	5 holes/left	94mm	
08.311.01.100033806	6 holes/left	108mm	
08.311.01.100033807	7 holes/left	123mm	
08.311.01.100033808	8 holes/left	135mm	
08.311.01.100044803	3 holes/right	69mm	
08.311.01.100044804	4 holes/right	81mm	
08.311.01.100044805	5 holes/right	94mm	
08.311.01.100044806	6 holes/right	108mm	
08.311.01.100044807	7 holes/right	123mm	
08.311.01.100044808	8 holes/right	135mm	

### Clavicle Locking Compression Plate-3



**Screws:**  $\Phi$ 2.4 Self- tapping locking screw- 3001824xx  
 $\Phi$ 2.7 Self- tapping locking screw- 3001827xx  
 $\Phi$ 3.5 Self- tapping locking screw- 1009835xx  
 $\Phi$ 3.5 Cortex screw- 61335xx  
 $\Phi$ 4.0 Cancellous screw- 61540xx

**Indication:** It is used for fracture on the apex of clavicle or combined Shaft fracture.



Ref #	Size	Length	Material
08.311.01.100033703	3 holes/left	69mm	T
08.311.01.100033704	4 holes/left	81mm	
08.311.01.100033705	5 holes/left	94mm	
08.311.01.100033706	6 holes/left	108mm	
08.311.01.100033707	7 holes/left	123mm	
08.311.01.100033708	8 holes/left	135mm	
08.311.01.100044703	3 holes/right	69mm	
08.311.01.100044704	4 holes/right	81mm	
08.311.01.100044705	5 holes/right	94mm	
08.311.01.100044706	6 holes/right	108mm	
08.311.01.100044707	7 holes/right	123mm	
08.311.01.100044708	8 holes/right	135mm	