




Solutions for Fibula Fracture





□ Distal Lateral Fibula Locking Compression Plate-6



Screws: Φ 2.4 Self- tapping locking screw- 3001824xx
 Φ 2.7 Self- tapping locking screw- 3001827xx
 Φ 3.5 Self- tapping locking screw- 1009835xx
 Φ 3.5 Cortex screw- 61335xx
 Φ 4.0 Cancellous screw- 61540xx

Indication: It is used for distal Fibula fracture.

Ref #	Size	Length	Material
08.311.01.101430803	3 holes/left	57.7mm	A
08.311.01.101430805	5 holes/left	77.7mm	
08.311.01.101430807	7 holes/left	97.7mm	
08.311.01.101430809	9 holes/left	117.7mm	
08.311.01.101430811	11 holes/left	137.7mm	
08.311.01.101431803	3 holes/right	57.7mm	
08.311.01.101431805	5 holes/right	77.7mm	
08.311.01.101431807	7 holes/right	97.7mm	
08.311.01.101431809	9 holes/right	117.7mm	
08.311.01.101431811	11 holes/right	137.7mm	

□ Distal Lateral Fibula Locking Compression Plate-3(2.7 locking screw on head)



Screws: Φ 2.4 Self- tapping locking screw- 3001824xx
 Φ 2.7 Self- tapping locking screw- 3001827xx
 Φ 3.5 Self- tapping locking screw- 1009835xx
 Φ 3.5 Cortex screw- 61335xx
 Φ 4.0 Cancellous screw- 61540xx

Indication: It is used for distal Fibula fracture.



Ref #	Size	Length	Material
08.311.01.100619803	3 holes/left	71.5mm	A
08.311.01.100619804	4 holes/left	84.5mm	
08.311.01.100619805	5 holes/left	96.5mm	
08.311.01.100619806	6 holes/left	109mm	
08.311.01.100619807	7 holes/left	121.5mm	
08.311.01.100618803	3 holes/right	71.5mm	
08.311.01.100618804	4 holes/right	84.5mm	
08.311.01.100618805	5 holes/right	96.5mm	
08.311.01.100618806	6 holes/right	109mm	
08.311.01.100618807	7 holes/right	121.5mm	

Distal Lateral Fibula Locking Compression Plate-2



Screws: Φ 2.4 Self-tapping locking screw- 3001824xx
 Φ 2.7 Self-tapping locking screw- 3001827xx
 Φ 3.5 Self-tapping locking screw- 1009835xx
 Φ 3.5 Cortex screw- 61335xx
 Φ 4.0 Cancellous screw- 61540xx

Indication: It is used for distal Fibula fracture and posterior lateral approach.

Ref #	Size	Length	Material
08.311.01.100611803	3 holes/left	77.5mm	A
08.311.01.100611804	4 holes/left	90.5mm	
08.311.01.100611805	5 holes/left	103.5mm	
08.311.01.100611806	6 holes/left	116.5mm	
08.311.01.100611807	7 holes/left	129.5mm	
08.311.01.100611808	8 holes/left	142.5mm	
08.311.01.100611809	9 holes/left	155.5mm	
08.311.01.100611810	10 holes/left	168.5mm	
08.311.01.100611811	11 holes/left	181.5mm	
08.311.01.100611812	12 holes/left	194.5mm	
08.311.01.100612803	3 holes/right	77.5mm	
08.311.01.100612804	4 holes/right	90.5mm	
08.311.01.100612805	5 holes/right	103.5mm	
08.311.01.100612806	6 holes/right	116.5mm	
08.311.01.100612807	7 holes/right	129.5mm	
08.311.01.100612808	8 holes/right	142.5mm	
08.311.01.100612809	9 holes/right	155.5mm	
08.311.01.100612810	10 holes/right	168.5mm	
08.311.01.100612811	11 holes/right	181.5mm	
08.311.01.100612812	12 holes/right	194.5mm	