

### Distal Anterior Lateral Tibia Locking Compression Plate



**Screws:**  $\Phi$  3.5 Self-tapping locking screw - 1009835xx  
 $\Phi$  3.5 Cortex screw- 61335xx  
 $\Phi$  4.0 Cancellous screw- 61540xx

**Indication:** It is used for distal tibia fracture.



| Ref #               | Size           | Length | Material |
|---------------------|----------------|--------|----------|
| 08.311.01.100361805 | 5 holes/left   | 80mm   | A        |
| 08.311.01.100361807 | 7 holes/left   | 106mm  |          |
| 08.311.01.100361809 | 9 holes/left   | 132mm  |          |
| 08.311.01.100361811 | 11 holes/left  | 158mm  |          |
| 08.311.01.100361813 | 13 holes/left  | 184mm  |          |
| 08.311.01.100361815 | 15 holes/left  | 210mm  |          |
| 08.311.01.100361817 | 17 holes/left  | 236mm  |          |
| 08.311.01.100361819 | 19 holes/left  | 262mm  |          |
| 08.311.01.100372805 | 5 holes/right  | 80mm   |          |
| 08.311.01.100372807 | 7 holes/right  | 106mm  |          |
| 08.311.01.100372809 | 9 holes/right  | 132mm  |          |
| 08.311.01.100372811 | 11 holes/right | 158mm  |          |
| 08.311.01.100372813 | 13 holes/right | 184mm  |          |
| 08.311.01.100372815 | 15 holes/right | 210mm  |          |
| 08.311.01.100372817 | 17 holes/right | 236mm  |          |
| 08.311.01.100372819 | 19 holes/right | 262mm  |          |

### Distal Medial Tibia Locking Compression Plate-1



**Screws:**  $\Phi$  3.5 Self-tapping locking screw - 1009835xx  
 $\Phi$  3.5 Cortex screw- 61335xx  
 $\Phi$  4.0 Cancellous screw- 61540xx

**Indication:** It is used for distal tibia fracture.



| Ref #               | Size           | Length | Material |
|---------------------|----------------|--------|----------|
| 08.311.01.100381803 | 3 holes/left   | 97mm   | A        |
| 08.311.01.100381805 | 5 holes/left   | 123mm  |          |
| 08.311.01.100381807 | 7 holes/left   | 149mm  |          |
| 08.311.01.100381809 | 9 holes/left   | 175mm  |          |
| 08.311.01.100381811 | 11 holes/left  | 201mm  |          |
| 08.311.01.100381813 | 13 holes/left  | 227mm  |          |
| 08.311.01.100381815 | 15 holes/left  | 253mm  |          |
| 08.311.01.100392803 | 3 holes/right  | 97mm   |          |
| 08.311.01.100392805 | 5 holes/right  | 123mm  |          |
| 08.311.01.100392807 | 7 holes/right  | 149mm  |          |
| 08.311.01.100392809 | 9 holes/right  | 175mm  |          |
| 08.311.01.100392811 | 11 holes/right | 201mm  |          |
| 08.311.01.100392813 | 13 holes/right | 227mm  |          |
| 08.311.01.100392815 | 15 holes/right | 253mm  |          |

### Q Distal Medial Tibia Locking Compression Plate-1



**Screws:**  $\Phi$ 3.5 Self- tapping locking screw- 1009835xx  
 $\Phi$ 3.5 Cortex screw- 61335xx  
 $\Phi$ 4.0 Cancellous screw- 61540xx

**Indication:** It is used for distal tibia fracture.

| Ref #               | Size           | Length | Material |
|---------------------|----------------|--------|----------|
| 08.311.01.100381703 | 3 holes/left   | 97mm   | T        |
| 08.311.01.100381705 | 5 holes/left   | 123mm  |          |
| 08.311.01.100381707 | 7 holes/left   | 149mm  |          |
| 08.311.01.100381709 | 9 holes/left   | 175mm  |          |
| 08.311.01.100381711 | 11 holes/left  | 201mm  |          |
| 08.311.01.100381713 | 13 holes/left  | 227mm  |          |
| 08.311.01.100381715 | 15 holes/left  | 253mm  |          |
| 08.311.01.100392703 | 3 holes/right  | 97mm   |          |
| 08.311.01.100392705 | 5 holes/right  | 123mm  |          |
| 08.311.01.100392707 | 7 holes/right  | 149mm  |          |
| 08.311.01.100392709 | 9 holes/right  | 175mm  |          |
| 08.311.01.100392711 | 11 holes/right | 201mm  |          |
| 08.311.01.100392713 | 13 holes/right | 227mm  |          |
| 08.311.01.100392715 | 15 holes/right | 253mm  |          |

### Q Locking T-Plates-1, oblique



**Screws:**  $\Phi$ 2.7 Self- tapping locking screw- 3001827xx  
 $\Phi$ 3.5 Self- tapping locking screw- 1009835xx

**Indication:** It is used for posterior ankle fracture.

| Ref #               | Size          | Length | Material |
|---------------------|---------------|--------|----------|
| 08.311.01.100519702 | 2 holes/left  | 57mm   | T        |
| 08.311.01.100519703 | 3 holes/left  | 70mm   |          |
| 08.311.01.100519704 | 4 holes/left  | 83mm   |          |
| 08.311.01.100519705 | 5 holes/left  | 96mm   |          |
| 08.311.01.100519706 | 6 holes/left  | 109mm  |          |
| 08.311.01.100528702 | 2 holes/right | 57mm   |          |
| 08.311.01.100528703 | 3 holes/right | 70mm   |          |
| 08.311.01.100528704 | 4 holes/right | 83mm   |          |
| 08.311.01.100528705 | 5 holes/right | 96mm   |          |
| 08.311.01.100528706 | 6 holes/right | 109mm  |          |