

Ø3.0mm Cannulated Screw(Countersunk head/Short thread)




With drill bit: cannulated, ϕ 2.0mm
(WG REF #30402103)

Indication: It is used for the fracture on metacarpal ,
Carpus and Plantar phalanges.

| Ref # | Size | Material |
|---------------------|----------|----------|
| 08.301.11.601133010 | 3.0x10mm | T |
| 08.301.11.601133012 | 3.0x12mm | |
| 08.301.11.601133014 | 3.0x14mm | |
| 08.301.11.601133016 | 3.0x16mm | |
| 08.301.11.601133018 | 3.0x18mm | |
| 08.301.11.601133020 | 3.0x20mm | |
| 08.301.11.601133022 | 3.0x22mm | |
| 08.301.11.601133024 | 3.0x24mm | |
| 08.301.11.601133026 | 3.0x26mm | |
| 08.301.11.601133028 | 3.0x28mm | |
| 08.301.11.601133030 | 3.0x30mm | |
| 08.301.11.601133032 | 3.0x32mm | |
| 08.301.11.601133034 | 3.0x34mm | |
| 08.301.11.601133036 | 3.0x36mm | |
| 08.301.11.601133038 | 3.0x38mm | |
| 08.301.11.601133040 | 3.0x40mm | |

Ø3.0mm Cannulated Screw(Countersunk head/half thread)




With drill bit: cannulated, ϕ 2.0mm
(WG REF #30402103)

Indication: It is used for the fracture on metacarpal ,
Carpus and Plantar phalanges.

| Ref # | Size | Material |
|---------------------|----------|----------|
| 08.301.11.601153010 | 3.0x10mm | T |
| 08.301.11.601153012 | 3.0x12mm | |
| 08.301.11.601153014 | 3.0x14mm | |
| 08.301.11.601153016 | 3.0x16mm | |
| 08.301.11.601153018 | 3.0x18mm | |
| 08.301.11.601153020 | 3.0x20mm | |
| 08.301.11.601153022 | 3.0x22mm | |
| 08.301.11.601153024 | 3.0x24mm | |
| 08.301.11.601153026 | 3.0x26mm | |
| 08.301.11.601153028 | 3.0x28mm | |
| 08.301.11.601153030 | 3.0x30mm | |
| 08.301.11.601153032 | 3.0x32mm | |
| 08.301.11.601153034 | 3.0x34mm | |
| 08.301.11.601153036 | 3.0x36mm | |
| 08.301.11.601153038 | 3.0x38mm | |
| 08.301.11.601153040 | 3.0x40mm | |



Ø4.5mm Cannulated Screw (Countersunk head)



With drill bit: cannulated, Φ 3.2mm
(WG REF #30402700)

Indication: It is used for the fracture on tibia plateau, anterior/posterior/lateral/medial ankle, proximal humerus.

| Ref # | Size | Material |
|---------------------|----------|----------|
| 08.301.11.601354520 | 4.5x20mm | T |
| 08.301.11.601354522 | 4.5x22mm | |
| 08.301.11.601354524 | 4.5x24mm | |
| 08.301.11.601354526 | 4.5x26mm | |
| 08.301.11.601354528 | 4.5x28mm | |
| 08.301.11.601354530 | 4.5x30mm | |
| 08.301.11.601354532 | 4.5x32mm | |
| 08.301.11.601354534 | 4.5x34mm | |
| 08.301.11.601354536 | 4.5x36mm | |
| 08.301.11.601354538 | 4.5x38mm | |
| 08.301.11.601354540 | 4.5x40mm | |
| 08.301.11.601354542 | 4.5x42mm | |
| 08.301.11.601354544 | 4.5x44mm | |
| 08.301.11.601354546 | 4.5x46mm | |
| 08.301.11.601354548 | 4.5x48mm | |
| 08.301.11.601354550 | 4.5x50mm | |
| 08.301.11.601354552 | 4.5x52mm | |
| 08.301.11.601354554 | 4.5x54mm | |
| 08.301.11.601354556 | 4.5x56mm | |
| 08.301.11.601354558 | 4.5x58mm | |
| 08.301.11.601354560 | 4.5x60mm | |
| 08.301.11.601354564 | 4.5x64mm | |
| 08.301.11.601354568 | 4.5x68mm | |
| 08.301.11.601354572 | 4.5x72mm | |



Ø3.5mm Cannulated Screw (Ball head)



With drill bit: cannulated, ϕ 2.7mm
(WG REF #30400712)

Indication: It is used for the fracture on tibia plateau, anterior/posterior/lateral/medial ankle, proximal humerus.

| Ref # | Size | Material |
|---------------------|----------|----------|
| 08.301.11.600233510 | 3.5x10mm | T |
| 08.301.11.600233512 | 3.5x12mm | |
| 08.301.11.600233514 | 3.5x14mm | |
| 08.301.11.600233516 | 3.5x16mm | |
| 08.301.11.600233518 | 3.5x18mm | |
| 08.301.11.600233520 | 3.5x20mm | |
| 08.301.11.600233522 | 3.5x22mm | |
| 08.301.11.600233524 | 3.5x24mm | |
| 08.301.11.600233526 | 3.5x26mm | |
| 08.301.11.600233528 | 3.5x28mm | |
| 08.301.11.600233530 | 3.5x30mm | |
| 08.301.11.600233532 | 3.5x32mm | |
| 08.301.11.600233534 | 3.5x34mm | |
| 08.301.11.600233536 | 3.5x36mm | |
| 08.301.11.600233538 | 3.5x38mm | |
| 08.301.11.600233540 | 3.5x40mm | |
| 08.301.11.600233542 | 3.5x42mm | |
| 08.301.11.600233544 | 3.5x44mm | |
| 08.301.11.600233546 | 3.5x46mm | |
| 08.301.11.600233548 | 3.5x48mm | |
| 08.301.11.600233550 | 3.5x50mm | |



Ø4.5mm Cannulated Screw (Ball head)



With drill bit: cannulated, ϕ 3.2mm
(WG REF #6090004)

Indication: It is used for the fracture on tibia plateau, anterior/posterior/lateral/medial ankle, proximal humerus.

| Ref # | Size | Material |
|---------------------|----------|----------|
| 08.301.01.600074020 | 4.5x20mm | T |
| 08.301.01.600074022 | 4.5x22mm | |
| 08.301.01.600074024 | 4.5x24mm | |
| 08.301.01.600074026 | 4.5x26mm | |
| 08.301.01.600074028 | 4.5x28mm | |
| 08.301.01.600074030 | 4.5x30mm | |
| 08.301.01.600074032 | 4.5x32mm | |
| 08.301.01.600074034 | 4.5x34mm | |
| 08.301.01.600074036 | 4.5x36mm | |
| 08.301.01.600074038 | 4.5x38mm | |
| 08.301.01.600074040 | 4.5x40mm | |
| 08.301.01.600074042 | 4.5x42mm | |
| 08.301.01.600074044 | 4.5x44mm | |
| 08.301.01.600074046 | 4.5x46mm | |
| 08.301.01.600074048 | 4.5x48mm | |
| 08.301.01.600074050 | 4.5x50mm | |
| 08.301.01.600074052 | 4.5x52mm | |
| 08.301.01.600074054 | 4.5x54mm | |
| 08.301.01.600074056 | 4.5x56mm | |
| 08.301.01.600074058 | 4.5x58mm | |
| 08.301.01.600074060 | 4.5x60mm | |
| 08.301.01.600074064 | 4.5x64mm | |
| 08.301.01.600074068 | 4.5x68mm | |
| 08.301.01.600074072 | 4.5x72mm | |

Ø7.3mm Cannulated Screw (Ball head/long thread)


With drill bit: cannulated, ϕ 5.0mm
(WG REF #6090202)

Indication: It is used for the fracture on femur neck and the fusion on tibia and astragalus and calcaneus.

| Ref # | l | Size | Material |
|---------------------|----|-----------|----------|
| 08.301.01.600058050 | 32 | 7.3x50mm | T |
| 08.301.01.600058055 | 32 | 7.3x55mm | |
| 08.301.01.600058060 | 32 | 7.3x60mm | |
| 08.301.01.600058065 | 32 | 7.3x65mm | |
| 08.301.01.600058070 | 32 | 7.3x70mm | |
| 08.301.01.600058075 | 32 | 7.3x75mm | |
| 08.301.01.600058080 | 32 | 7.3x80mm | |
| 08.301.01.600058085 | 32 | 7.3x85mm | |
| 08.301.01.600058090 | 32 | 7.3x90mm | |
| 08.301.01.600058095 | 32 | 7.3x95mm | |
| 08.301.01.600058100 | 32 | 7.3x100mm | |
| 08.301.01.600058105 | 32 | 7.3x105mm | |
| 08.301.01.600058110 | 32 | 7.3x110mm | |
| 08.301.01.600058115 | 32 | 7.3x115mm | |
| 08.301.01.600058120 | 32 | 7.3x120mm | |

Ø7.3mm Cannulated Screw (Ball head/short thread)


With drill bit: cannulated, ϕ 5.0mm
(WG REF #6090202)

Indication: It is used for the fracture on femur neck and the fusion on tibia and astragalus and calcaneus.

| Ref # | l | Size | Material |
|---------------------|----|-----------|----------|
| 08.301.01.600048050 | 16 | 7.3x50mm | T |
| 08.301.01.600048055 | 16 | 7.3x55mm | |
| 08.301.01.600048060 | 16 | 7.3x60mm | |
| 08.301.01.600048065 | 16 | 7.3x65mm | |
| 08.301.01.600048070 | 16 | 7.3x70mm | |
| 08.301.01.600048075 | 16 | 7.3x75mm | |
| 08.301.01.600048080 | 16 | 7.3x80mm | |
| 08.301.01.600048085 | 16 | 7.3x85mm | |
| 08.301.01.600048090 | 16 | 7.3x90mm | |
| 08.301.01.600048095 | 16 | 7.3x95mm | |
| 08.301.01.600048100 | 16 | 7.3x100mm | |
| 08.301.01.600048105 | 16 | 7.3x105mm | |
| 08.301.01.600048110 | 16 | 7.3x110mm | |
| 08.301.01.600048115 | 16 | 7.3x115mm | |
| 08.301.01.600048120 | 16 | 7.3x120mm | |