



Solutions for Cannulated Screws





Ø3.0mm Self- tapping Cannulated Screw



With drill bit: Φ 2.0 Cannulated drill (WG REF #30402206)
 Φ 2.0 conical drill (WG REF# 30402506)

Indication: It is used for the fracture on metacarpal ,
 Carpus and Plantar phalanges.

Ref #	Size	Material
08.301.11.601573008	3.0x8mm	T
08.301.11.601573010	3.0x10mm	
08.301.11.601573012	3.0x12mm	
08.301.11.601573014	3.0x14mm	
08.301.11.601573016	3.0x16mm	
08.301.11.601573018	3.0x18mm	
08.301.11.601573020	3.0x20mm	

Ø3.5mm Self- tapping Cannulated Screw




With drill bit: Φ 2.5 cannulated bone drill (WG REF #30402106)
 Φ 2.5 Conical bone drill (WG REF #30402406)

Indication: It is used for the fracture on metacarpal ,
 Carpus and Plantar phalanges.

Ref #	Size	Material
08.301.11.601613516	3.5x16mm	T
08.301.11.601613518	3.5x18mm	
08.301.11.601613520	3.5x20mm	
08.301.11.601613522	3.5x22mm	
08.301.11.601613524	3.5x24mm	
08.301.11.601613526	3.5x26mm	
08.301.11.601613528	3.5x28mm	
08.301.11.601613530	3.5x30mm	

Ø4.0mm Self- tapping Cannulated Screw




With drill bit: Φ 3.0 canulated bone drill (WG REF #30402006)
 Φ 3.0 Conical bone drill (WG REF #30402306)

Indication: It is used for the fracture on metacarpal , Carpus and Plantar phalanges.

Ref #	Size	Material
08.301.11.601654016	4.0x16mm	T
08.301.11.601654018	4.0x18mm	
08.301.11.601654020	4.0x20mm	
08.301.11.601654022	4.0x22mm	
08.301.11.601654024	4.0x24mm	
08.301.11.601654026	4.0x26mm	
08.301.11.601654028	4.0x28mm	
08.301.11.601654030	4.0x30mm	

Ø4.5mm Self- tapping Cannulated Screw




With drill bit: Φ 4.6- 3.2 drill (WG REF #30402906)
 Φ 4.4- 3 drill (WG REF #30403006)

Indication: It is used for the fracture on tibia plateau, anterior/posterior/lateral/medial ankle, proximal humerus.

Ref #	Size	Material
08.301.11.601694524	4.5x24mm	T
08.301.11.601694526	4.5x26mm	
08.301.11.601694528	4.5x28mm	
08.301.11.601694530	4.5x30mm	
08.301.11.601694532	4.5x32mm	
08.301.11.601694534	4.5x34mm	
08.301.11.601694536	4.5x36mm	
08.301.11.601694538	4.5x38mm	
08.301.11.601694540	4.5x40mm	
08.301.11.601694542	4.5x42mm	
08.301.11.601694544	4.5x44mm	
08.301.11.601694546	4.5x46mm	
08.301.11.601694548	4.5x48mm	
08.301.11.601694550	4.5x50mm	



Ø5.0mm Self- tapping Cannulated Screw




With drill bit: Φ 4-5.6 drill (WG REF #30403606)
 Φ 4.5- 6 drill (WG REF # 30403706)

Indication: It is used for the fracture on tibia plateau, anterior/posterior/lateral/medial ankle, proximal humerus.

Ref #	Size	Material
08.301.11.601735035	5.0x35mm	T
08.301.11.601735040	5.0x40mm	
08.301.11.601735045	5.0x45mm	
08.301.11.601735050	5.0x50mm	
08.301.11.601735055	5.0x55mm	
08.301.11.601735060	5.0x60mm	
08.301.11.601735065	5.0x65mm	
08.301.11.601735070	5.0x70mm	
08.301.11.601735075	5.0x75mm	
08.301.11.601735080	5.0x80mm	



Ø7.0mm Self-tapping Cannulated Screw



With drill bit: Φ 4.5- 6.7
 Reamer 1 (WG REF#30404406)
 Reamer 2 (WG REF#30404506)
 Reamer 3 (WG REF#30404606)
 Φ 5.0- 6.7
 Reamer 1 (WG REF# 30404706)
 Reamer 2 (WG REF# 30404806)
 Reamer 3 (WG REF# 30404906)

Indication: It is used for the fusion of tibia and astragalus.

Ref #	Size	Material
08.301.11.601777040	7.0x40mm	T
08.301.11.601777045	7.0x45mm	
08.301.11.601777050	7.0x50mm	
08.301.11.601777055	7.0x55mm	
08.301.11.601777060	7.0x60mm	
08.301.11.601777065	7.0x65mm	
08.301.11.601777070	7.0x70mm	
08.301.11.601777075	7.0x75mm	
08.301.11.601777080	7.0x80mm	
08.301.11.601777085	7.0x85mm	
08.301.11.601777090	7.0x90mm	
08.301.11.601777095	7.0x95mm	
08.301.11.601777100	7.0x100mm	